



easy-graft[®] **CLASSIC**
synthetic bone grafting system



More Consistent, Predictable Results

If the ability to place an implant is your primary measure of success, GUIDOR® *easy-graft*® provides a predictable way to maintain alveolar bone for rehabilitation with dental implants after tooth extraction.⁵

Clinical Case: Rigid, Porous Scaffold of GUIDOR *easy-graft* Achieves Implant Success

A 23-year-old female patient without medical contraindication for implant therapy presented with a maxillary right second premolar that due to extensive caries could not be saved. The treatment plan included a delayed implant placement protocol.



Figure 1: Pre-op Clinical View



Figure 2: Pre-op Radiograph



Figure 3: Atraumatic Extraction



Figure 4: Socket Filled with GUIDOR *easy-graft*



Figure 5: Clinical Result 3 Weeks Post-op
Gradual proliferation of the epithelium over the grafted site.



Figure 6: Clinical Result 4 Months Post-Op
Grafted area is covered with newly formed keratinized epithelium.



Figure 7: Newly Formed Bone at 4 Month Re-entry
Adequate preservation of ridge dimension can be seen.




Figure 8: Histological Analysis
*GUIDOR *easy-graft* particles surrounded by newly formed bone.*
Key: EG = GUIDOR *easy-graft*
NB = newly formed bone



Figure 9: Implant Placement



Figure 10: Clinical Result After 4 Months Post-implant Placement



Figure 11: Clinical View After Crown Placement



Figure 12: Final Radiograph at 16 Months
No residual particles of the grafting material are detected and the implant seems surrounded by vital bone.




easy-graft® CLASSIC synthetic bone grafting system

GUIDOR® easy-graft® is a patented¹ bone grafting solution offering a fully resorbable, synthetic bone graft substitute to clinicians who seek stability, osteoconductivity and control of placement. With over 250,000 applications worldwide², GUIDOR *easy-graft* has established itself as a leading synthetic bone grafting material for your ridge preservation and other bone regenerative needs.


Designed to Make It Easy

GUIDOR *easy-graft* is a complete system that helps clinicians provide more predictable clinical results. Each system is comprised of β -tricalcium phosphate (β -TCP) granules coated with a poly(lactide-co-glycolide) (PLGA) that are mixed with N-methyl-2-pyrrolidone (NMP) liquid activator called BioLinker® to form a permeable, moldable material which hardens to form a stable, porous scaffold.


This unique combination of materials in GUIDOR *easy-graft* provides the foundation for easy preparation, dispensing and shaping. The porosity of the granule supports material resorption and bone regeneration.




The system contains BioLinker® and a β -TCP granule with polymer coating.



When mixed, BioLinker softens the polymer coating creating a sticky surface that can be compressed and shaped.



BioLinker is flushed out of the material when in contact with body fluids.



GUIDOR *easy-graft* hardens in approximately one minute into a stable, porous scaffold of interconnected granules.

Soft from the Syringe. Hard in the Defect.



Placing bone graft material doesn't have to be complex. In fact, a majority of clinicians surveyed said they use GUIDOR *easy-graft* because preparation and placement are simpler and faster than other bone graft materials² and because use of GUIDOR *easy-graft* does not require a membrane in many cases.³

- Easy-to-handle; syringe directly into the defect
- Material hardens into a stable, porous scaffold in approximately one minute
- Ideal for ridge preservation and filling voids around immediate implant placements

1. US Patent: 7,731,756; 8,153,148; 8,163,030

2. Data on file. Degradable Solutions AG. 3. GUIDOR *easy-graft* Customer Satisfaction Survey, February 2014.

Simple Steps for Minimally Invasive Ridge Preservation


Prep.



1
Atraumatic Extraction



2 KEY STEP
Socket Preparation
Thoroughly clean alveolar walls.
Cause bleeding.



3
Socket Evaluation



4 KEY STEP
Create Additional Retention
In mesial and distal walls.

Dispense.



5
Prepare BioLinker®



6
Introduce BioLinker



7
Wet Granules with BioLinker



8
Expel Excess BioLinker


Shape.



9
Inject Material into the Site




10 KEY STEP
Compress Material Firmly



11
Shape and Contour
GUIDOR® *easy-graft*® hardens in contact with body fluids in minutes.

Clinical Preferences⁴ when using GUIDOR® *easy-graft*® CLASSIC synthetic bone grafting system



- Prior to placement of GUIDOR® *easy-graft*® CLASSIC, **95%** of respondents **cleaned the defect with mechanical devices**, such as curette and burs.
- **52%** of respondents did **not cover the site** with anything after placing GUIDOR® *easy-graft*® CLASSIC in the defect.
- **83%** of respondents **placed sutures** over the site.

Opinions after use⁴


An independent healthcare market research firm conducted a blind product evaluation, where 80 US dentists tried GUIDOR® *easy-graft*® CLASSIC Synthetic Bone Grafting System for the first time. From 2014-2015, the respondents indicated they completed a total of 300+ cases with an average of 5 cases per clinician and provided feedback after 3 months (N=75) and 6 months (N=65).⁴

UP TO 90%

respondents were satisfied with the clinical outcomes when using GUIDOR® *easy-graft*® CLASSIC


Satisfaction with Clinical Outcomes

% of Clinicians



Agreement with Positive Statements

% of Clinicians



90%

confirmed that GUIDOR® *easy-graft*® CLASSIC simplifies bone grafting procedures

95%

of the respondents did not report any patient complaints after using GUIDOR® *easy-graft*® CLASSIC

Any complications?






97%

of the respondents did not report any signs of infection or inflammatory reactions after GUIDOR® *easy-graft*® CLASSIC

GUIDOR® *easy-graft*® Synthetic Bone Grafting System


GUIDOR *easy-graft* is a fully resorbable bone grafting system that contains a syringe filled with coated granules and an ampule of the liquid activator, BioLinker®. Once mixed, GUIDOR *easy-graft* can be dispensed and shaped directly into the defect where it will harden into a stable, porous scaffold in approximately one minute.

Product Size	Item #	Granule Size	Content	Color Code
Large	C11-008	500-1000µm	3 systems x 0.4ml	
Medium	C11-078	500-1000µm	3 systems x 0.25ml	
Small	C11-018	500-630µm	3 systems x 0.15ml	

1 system = 1 syringe of GUIDOR *easy-graft* granules and 1 ampule of BioLinker®

Size Selector: The Amount of GUIDOR *easy-graft* needed to fill a socket

The following estimates are based on GUIDOR *easy-graft* placement in dentiform sockets. Material needs in clinical use may vary.



Indications for Use

GUIDOR *easy-graft* is indicated for the treatment of intraoral/maxillofacial osseous defects.

Dental and maxillofacial indications may include:

- Extraction defects (alveolar ridge preservation)
- Periodontal defects
- Peri-implant defects
- Augmentation of deficient alveolar crest (e.g. Guided Bone Regeneration)
- Sinus floor augmentation
- Defects after surgical extractions
- Defects after removal of bone cysts
- Defects after root resection of apicoectomy
- Defects after removal of autologous bone

Contraindications

GUIDOR *easy-graft* should not be used in pregnant or nursing women.

Possible Adverse Effects

Possible adverse reactions associated with the use of the device include: eye, respiratory and skin irritation.

